TECHNICAL DATA SHEET



Antiquax Original Wax Polish -Aerosol

Waterloo Road, London NW2 7TX Telephone: +44(0)20 8450 4666

WA

Туре	A very high quality traditional type polish based upon premium performance waxes, carnauba and beeswax, reinforced with a blend of hydrocarbon waxes to aid application and improve durability. This product has been specially modified to deliver Antiquax Original Wax Polish in an easy to use form.
Intended Use	For regular use on all fine natural wood and antique furniture.
Special Features	Improved ease of application. Provides a long lasting durable wax finish with soft mellow sheen to enhance grain and patina. Silicone free. Resists dust and fingerprinting.
Preparation	The surface must be clean and dry. Antiquax Original Wax Polish can be used on previously waxed surfaces. For best results initially clean the surface with Antiquax Furniture Cleaner .
Application	Shake thoroughly before use. Spray apply sparingly and work into the surface with the grain using a soft lint free cloth. Work well into the grain and allow to dry to a light haze. Remove by polishing in a circular motion using a clean dry cloth.
Coats	One coat is usually sufficient but a second coat may be applied to obtain a higher sheen.
Drying Time	Quick drying.
Re-coating Time	Re-apply if required after a few minutes.
Application Temp.	Ideally 15-25°C
Coverage	N/A
Maintenance	Re-apply when necessary.
Compatibility	Compatible with most timber finishes. If in doubt check on an inconspicuous test area. For best results and to prevent a blotchy appearance, timber should be previously sealed with a varnish, oil or sanding sealer.
Settlement	Shake thoroughly before use
Shelf-life/Storage	A very long time if stored in original sealed container in a cool dry place.
Hazards	Extremely flammable. Irritating to eyes. Vapours may cause drowsiness. Keep away from children. Keep away from all sources of ignition. No smoking. A full Health & Safety Data Sheet is available on request or from www.antiquax.co.uk
Clean-Up	Hands and application cloths can be cleaned in warn soapy water.
Sizes	250 ml aerosols