



www.antiquax.co.uk

Waterloo Road, London NW2 7TX
Telephone: +44(0)20 8450 4666

TECHNICAL DATA SHEET

ANTIQUAX PRODUCT NAME
Silver Dip

CODE ANTQSILD250

Type	Acid free Solution
Intended Use	Antiquax Silver Dip is acid free formation that will remove tarnish from silver cutlery, jewellery and small silver items.
Special Features	No rubbing required. Cleans and shines in minutes. Suitable for engraved and filigree pieces.
Application	Immerse silver items in the jar until tarnish disappears (but not too longer than 5 mins.) Rinse with water and dry. If tarnish persists, re-dip and polish using an Antiquax Polishing Cloth.
Coats	N/A
Drying Time	N/A
Re-coating Time	N/A
Application Temp.	Warm water temperature.
Coverage	N/A
Maintenance	Re-treat when necessary
Compatibility	Antiquax is not suitable for use on iron, steel, stainless steel, nickel, worn silver or artificially oxidised silver. If in doubt test on a small area before proceeding with a full dip.
Settlement	N/A
Shelf-life/Storage	Store in a cool dry place.
Hazards	Can cause skin irritation. Can cause serious eye irritation. Suspected of causing Cancer. For details read the label carefully.
Clean-Up	Warm soapy water.
Sizes	250ml/

This information is provided in good faith based upon data and sources believed to be reliable and correct. Conditions of use, outside the control of the Company, dictate that no responsibility can be assumed for the ultimate performance of this product.